



UCLA OLYMPIC ANALYTICAL LABORATORY
DEPARTMENT OF PHARMACOLOGY
UCLA SCHOOL OF MEDICINE
2122 ORANVILLE AVENUE
LOS ANGELES, CALIFORNIA 90025
PHONE (310) 825-2635
FAX (310) 206-9077

CONFIDENTIAL

May 27, 2004

Terry Madden
The United States Anti-Doping Agency
2550 Tenderfoot Hill St., Suite 200
Colorado Springs, CO 80906-7346

RE: Specimen number [REDACTED]


Dear Mr. Madden

Please find enclosed the documentation package for the screen, the A confirmation, and the B confirmation on the case identified above.

Enclosed are computer generated images and graphs supporting our conclusion and the drug testing report.

Please feel free to call if you have any questions

Sincerely


Don H. Catlin, M.D.
Director

ASB 05/01/04



CONFIDENTIAL
DOCUMENTATION

SAMPLE IDENTIFICATION:

Organization requesting test: USADA
Date of sample collection: Mar 16, 2004
Site ID: OOC



Substance identified: recombinant human Erythropoietin (rHuEPO)

TABLE OF CONTENTS

	Page
CHAIN OF CUSTODY DOCUMENTATION FOR SAMPLE CONTAINER	
Official record, lab copy	4
Courier documentation	5
Box chain of custody documentation	6
Laboratory confirmation of receipt of sample	8
Individual seals chain of custody documentation	9
"A" Screen aliquot chain of custody documentation	10
"A" bottle chain of custody documentation	12
"A" confirmation aliquot chain of custody documentation	13
"A" sample stability test aliquot chain of custody documentation	15
"B" bottle chain of custody documentation	16
"B" confirmation aliquot chain of custody documentation	17
"B" sample stability test aliquot chain of custody documentation	18
Description of the EPO assay	19
"A" SAMPLE SCREENING DOCUMENTATION	
Figure 1: Electropherogram with isoform patterns - Includes this sample, Standards, Negative QC, Positive QC, and one other sample in batch	24
Results of Screen	25

"A" SAMPLE CONFIRMATION DOCUMENTATION

Figure 2 : Electropherogram with isoform patterns - Includes this sample, Standards, Negative QC, Positive QC 26

Figure 3 : Isoform profile A confirmation 27

Results of A confirmation 28

"B" SAMPLE CONFIRMATION DOCUMENTATION

Figure 4 : Electropherogram with isoform patterns - Includes this sample, Standards, Negative QC, Positive QC 29

Figure 5 : Isoform profile of B confirmation 30

Results of B confirmation 31

SAMPLE STABILITY DOCUMENTATION

Figure 6 : Electropherogram of stability test (A confirmation) 32

Figure 7 : Electropherogram of stability test (B confirmation) 32

Figure 8 : Overlay of isoform profiles of screen, A-, and B-confirmation 33



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

CONFIDENTIAL

UCLA OLYMPIC ANALYTICAL LABORATORY
Don S. Catlin, M.D., Director
2122 Granville Avenue, Los Angeles, CA 90025
PH (310) 825-2635 FAX (310) 206-9077

BATCH CHAIN OF CUSTODY

SHIPMENT RELEASED BY	TO	DATE/TIME	PURPOSE
<i>SPS</i> <i>A. H. Montoya</i> SIGNATURE SPINT - COURIER	<i>[Signature]</i> Fereshteh Geshad FILL OUT SECTION RIGHT BELOW	3/18/04 10:10AM	TRANSFER

Said shipment consists of:
Airbill # or Package Tracking # or description

- Bags
No. _____
(color) _____
- Boxes
No. _____
- Envelopes
No. _____
121A1W422 10080795 (A)
121A1W422 10083069 (A)
121A1W422 10132041 (A)
121A1W422 10133223 (A)

(other) explain: _____

INITIAL NEXT TO THE AIRBILL # TO INDICATE THE BATCH.


Integrity of the shipment is (check one):

undamaged damaged (describe: _____)

see reverse side for continuation.



UCLA Olympic Analytical Laboratory
 UCLA School of Medicine
 2122 Granville Ave Los Angeles CA 90025
 Phone (310) 825-2635 FAX (310) 206-9077



After inspection, custody of the shipment is as follows (document transfers of the shipment to Receiving Chemist, Room 118, Coldroom, etc.):

RELEASED BY	TO	DATE	PURPOSE
<i>[Signature]</i> Parsifal Delgado	REFRIGERATOR # 17	MAR 18 2004	STORAGE
REFRIGERATOR # 17	<i>[Signature]</i>	MAR 18 2004	Accessioning

The shipment is unpacked in the Accessioning Room and assigned a UCLA folder label:

Client: USADA Folder # [REDACTED] UCLA CODE: [REDACTED]
 The number of A samples is: [REDACTED] number of B samples is: [REDACTED] / FDN/18/04
 The number of Athlete Signature Forms/ Official Records is: 24
 The number of Specimen Custody Documents (2624) is: _____

After unpacking and batching the shipment, BOTTLE custody is as follows:

RELEASED BY	TO	DATE	PURPOSE
<i>[Signature]</i>			Aliquot A's in Rm 118
	COLDROOM	MAR 18 2004	STORAGE
			<u>GC 3/23/04</u> <u>51 31-100</u>

2150P\LC034PRM\CUSBRATR.2003

Form approved AR 2/2/03



UCLA Olympic Analytical Laboratory
 UCLA School of Medicine
 2122 Granville Ave Los Angeles CA 90025
 Phone (310) 825-2635 FAX (310) 206-9077



UNITED STATES ANTI-DOPING AGENCY
Sample Manifest

USADA 6012

Sample Information

- 1. DATE(S) OF COLLECTION: 03-16-04
(WRITE ALL DATES IN FULL LETTERS)
- 2. NAME OF EVENT (If OOC write "OOO"): OOO
- 3. SITE ID Number (If OOC write "OOO"): OOO
- 4. CITY & STATE where samples shipped from: Albuquerque, N.M.

5. TEST KIT NUMBERS SHIPPED:

[REDACTED]				

Lab Codes: no problems A=problem w/ bottle B=problem w/ bottle AB=problem with both bottles

- 6. COURIER: World Courier UPS Other
- 7. WAYBILL/Tracking NUMBER: 1Z 1A1 7V4 22 1013 327 3
- 8. BAG SEAL NUMBER (If OOC write "OOO"): OOO
- 9. DATE SEALED: 03 17 04
Month Day Year
- 10. TIME SEALED: 9:00 AM PM
- 11. Replacement Seal Number (if needed): _____
- 12. Reason for Replacement: _____
- 13. Date of Replacement: _____
Month Day Year
- 14. Time of Replacement: _____ AM PM
- 15. LEAD OOC: Mickey Williams Mickey Williams
Signature Official Name

For Laboratory Use Only

- LABORATORY CONFIRMATION OF RECEIPT OF SAMPLES: PLEASE COMPLETE & FAX TO USADA AT (719) 755-2001
- 1. DATE RECEIVED: 3 18 04
Month Day Year
 - 2. BAG SEAL NUMBER: 382
(If OOC write "OOO")
 - 3. COMMENTS: _____
 - 4. LABORATORY REPRESENTATIVE: [Signature] Fereshteh Deishad
Signature Printed Name



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

UCLA OLYMPIC ANALYTICAL LABORATORY

SAMPLE CUSTODY - DOCUMENTATION OF SEAL INTEGRITY

FOLDER: USADA [redacted]

UCLA CODE	BOTTLE NUMBER	A	B	* If not intact, describe below
[redacted]	[redacted]	0	0	
[redacted]	[redacted]	0	0	
[redacted]	[redacted]	X	X	
[redacted]	[redacted]	X	X	
[redacted]	[redacted]	X	X	
[redacted]	[redacted]	X	X	
[redacted]	[redacted]	X	X	
[redacted]	[redacted]	X	X	
[redacted]	0	0	0	
[redacted]	0	0	0	
[redacted]	0	0	0	

* Instructions: Verify by "X" that the custody of each specimen is intact. If not intact, record "NO" and describe in the space provided.
Note: Entries for internal QC are listed as "0".

The above samples were received in the condition stated above.

Signature [Signature] Date 3/18/04
Print Fereshteh Dalshad Time 12:50 PM
Verified (initials and date) ASJ 3/18/04

gc 3/18/04

SS 4/5/04
AB 4/5/04



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

CHAIN OF CUSTODY DOCUMENTATION

Circle: BOTTLE OR ALIQUOT(S) ACR B
 Aliquot: SCREEN (BATCH OR INDIVIDUAL) OR INDIVIDUAL CONFIRMATION
 Organization: USADA UCLA folder no. USADA
 If batch, UCLA code numbers: NA
 If individual, bottle no. [redacted] = UCLA code no. [redacted]
 If multiple replicate aliquots indicate here N = —
 Circle assay: I E CF FF DB TE ONTRAK DIU OTHER EPO
AC-BASE

IT IS REQUIRED TO PRINT EACH NAME AT LEAST ONCE ON THIS PAGE

Released by:	To	Date	Purpose
<i>Fereshteh Delshad</i> Fereshteh Delshad	REFRIGERATOR #22	MAR 1 8 2004	STORAGE
REFRIGERATOR #22	<i>A Range</i> Annie Ramsay	MAR 1 9 2004	Begin Assay
<i>A Range</i>	FREEZER #24	MAR 1 9 2004	STORAGE
FREEZER #24	<i>A Range</i>	MAR 2 4 2004	Continue Assay
<i>A Range</i>	Electrophoresis unit #2	MAR 2 4 2004	IEF
Electrophoresis unit #2	<i>A Range</i>	MAR 2 4 2004	Continue Assay
<i>A Range</i>	<i>Henry Truong</i> Henry Truong	MAR 2 4 2004	Continue Assay
<i>Henry Truong</i>	REFRIGERATOR #22	MAR 2 4 2004	Incubation
REFRIGERATOR #22	<i>A Range</i>	MAR 2 5 2004	Continue Assay
<i>A Range</i>	FLUORESCENCE CAMERA	MAR 2 5 2004	DATA ACQUISITION S. 415104

CHSOP/LOQSA/FRM-CUSTODY

Form approved *Kay - 1/13/04*



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

IT IS REQUIRED TO PRINT EACH NAME AT LEAST ONCE ON THIS PAGE

Released by:	To:	Date	Purpose
ELVOCHEM CAMERA	<i>J. Raup</i> Annie Ramsayer	MAR 2 5 2004	Transfer
<i>J. Raup</i>	long term STORAGE	MAR 2 5 2004	End Assay SD- 415104
			<i>J. Raup</i> 9/2/04

URSOP/LOGS&FRM/CLSTORY

Form approved

J. Raup 11/4/04
42 1/15/04



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

CHAIN OF CUSTODY DOCUMENTATION

Circle: BOTTLE OR ALIQUOT(S) A OR B

Aliquot: SCREEN (BATCH OR INDIVIDUAL) OR INDIVIDUAL CONFIRMATION

Organization: USADA UCLA folder no. USADA [REDACTED]

If batch, UCLA code numbers: NA

If individual, bottle no. [REDACTED] = UCLA code no. [REDACTED]

If multiple replicate aliquots indicate here N = —

Circle assay: I B CF FF BB TE ONTRAK DIU OTHER _____
AC-BASE

IT IS REQUIRED TO PRINT EACH NAME AT LEAST ONCE ON THIS PAGE

Released by:	To:	Date	Purpose
COLDROOM	Ferganah Delshad <i>[Signature]</i>	MAR 25 2004	Aliquot in Room 108
<i>[Signature]</i>	Positive Freezer <i>[Signature]</i>	MAR 25 2004	STORAGE S 415/04 AS 415/04

UCLA OLYMPIC ANALYTICAL LABORATORY

Form approved *[Signature]* - 4/2/04



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

CHAIN OF CUSTODY DOCUMENTATION

Circle: BOTTLE OR ALIQUOT(S) ADR B

Aliquot: SCREEN (BATCH OR INDIVIDUAL) OR INDIVIDUAL CONFIRMATION

Organization: USADA UCLA folder no. USADA

If batch, UCLA code numbers: NA

If individual, bottle no. [redacted] - UCLA code no. [redacted]

If multiple replicate aliquots indicate here N =

Circle assay: I II CF FF BB TE ONTRAK DIU OTHER EPD

IT IS REQUIRED TO PRINT EACH NAME AT LEAST ONCE ON THIS PAGE

Released by:	To:	Date	Purpose
<i>[Signature]</i> Fereshteh Delehdad	REFRIGERATOR # 22	MAR 2 5 2004	STORAGE
REFRIGERATOR #22	<i>[Signature]</i> Andreas Dreidbach	MAR 2 5 2004	Transfer
<i>[Signature]</i>	FREEZER #24	MAR 2 5 2004	STORAGE
FREEZER #24	<i>[Signature]</i> Henry Truong	MAR 2 9 2004	Begin Assay
<i>[Signature]</i>	REFRIGERATOR 22	MAR 2 9 2004	STORAGE
REFRIGERATOR #22	<i>[Signature]</i> Anne Bamseyer	MAR 3 0 2004	Continue Assay
<i>[Signature]</i>	Electrophoresis unit # 2	MAR 3 0 2004	IEF
Electrophoresis unit # 2	<i>[Signature]</i>	MAR 3 0 2004	Continue Assay
<i>[Signature]</i>	REFRIGERATOR #22	MAR 3 0 2004	Incubation
REFRIGERATOR #22	<i>[Signature]</i>	MAR 3 1 2004	Continue Assay 53-715704

G:\SPELOS&FRM\CHEMISTRY

Form approved *[Signature]* 1/19/04



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

IT IS REQUIRED TO PRINT EACH NAME AT LEAST ONCE ON THIS PAGE

Released by:	To:	Date	Purpose
<i>A. Ranga</i> Annie Ramaeyer	FLUOCHEM CAMERA	MAR 31 2004	DATA ACQUISITION
FLUOCHEM CAMERA	<i>A. Ranga</i>	MAR 31 2004	Transfer
<i>A. Ranga</i>	Long Term STORAGE	MAR 31 2004	End Assay 5.15.04
			<i>AS 4/2/04</i>

G:\SOP\LOGS\FRM\ICUSTODY

Form approved *[Signature]* 1/12/04



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

CHAIN OF CUSTODY DOCUMENTATION

Circle: BOTTLE OR ALIQUOT(S) A OR B

Aliquot: SCREEN (BATCH OR INDIVIDUAL) OR INDIVIDUAL CONFIRMATION

Organization: USADA UCLA folder no USADA [redacted]

If batch, UCLA code numbers: NA

If individual, bottle no. [redacted] = UCLA code no. [redacted]

If multiple replicate aliquots indicate here N = NA

Circle assay: I II CP FF BB TE ONTRAK DLI OTHER EPD Activity

IT IS REQUIRED TO PRINT EACH NAME AT LEAST ONCE ON THIS PAGE

Released by	To:	Date	Purpose
FREEZER #4	Henry Truong Henry Truong	MAR 29 2004	Aliquoting
Henry Truong	Water bath Room #125	MAR 29 2004	Incubation
WATER BATH ROOM #125	J Range Annie Flameeyer	MAR 30 2004	Continue Assay
J Range	unit #2 Electrophoresis	MAR 30 2004	IEF
Electrophoresis unit #2	J Range	MAR 30 2004	Continue Assay
J Range	REFRIGERATOR #22	MAR 30 2004	Incubation
REFRIGERATOR #22	J Range	MAR 31 2004	Continue Assay
J Range	FLUOCHEM CAMERA	MAR 31 2004	DATA ACQUISITION
FLUOCHEM CAMERA	J Range	MAR 31 2004	Transfer
J Range	long term STORAGE	MAR 31 2004	End Assay 5/15/04

0150F-1005&FORM CUSTODY

Form approved [Signature] 4/15/04

AS 4/2/04



UCLA Olympic Analytical Laboratory
 UCLA School of Medicine
 2122 Granville Ave Los Angeles CA 90025
 Phone (310) 825-2635 FAX (310) 208-9077



CHAIN OF CUSTODY DOCUMENTATION

Circled: BOTTLE OR ALIQUOT(S) A OR B

Aliquot: SCREEN (BATCH OR INDIVIDUAL) OR INDIVIDUAL CONFIRMATION

Organization: USADA UCLA folder no. USADA

If batch, UCLA code numbers: NA

If individual, bottle no. [REDACTED] = UCLA code no. [REDACTED]

If multiple replicate aliquots indicate here N =

Circled assay: I II CF PF BB TE ONTRAK DUJ OTHER
AC-BASE

IT IS REQUIRED TO PRINT EACH NAME AT LEAST ONCE ON THIS PAGE

Released by:	To:	Date	Purpose
	<i>[Signature]</i> Fereshteh Dolshad	MAR 25/2004	Transfer
<i>[Signature]</i>	Positive Freezer #2	MAR 25/2004	STORAGE
Positive Freezer #1	<i>[Signature]</i>	MAY 24 2004	Aliquot in Room 108
<i>[Signature]</i>	Positive Freezer #2	MAY 24 2004	STORAGE 5/25/04 F. S. 05/24/04

Q:\SDP\LOGS&FORM\CUSTODY

Form approved SD 7/15/02



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

CHAIN OF CUSTODY DOCUMENTATION

Circle BOTTLE OR ALIQUOT(S) A OR B
Aliquot: SCREEN (BATCH OR INDIVIDUAL) OR INDIVIDUAL CONFIRMATION
Organization: USADA UCLA folder no. USADA [redacted]
If batch, UCLA code numbers 2A
If individual, bottle no. [redacted] = UCLA code no. [redacted]
If multiple replicate aliquots indicate here N = 2A
Circle assay: I II CF FF FB TE ONTRAK DIU OTHER 1A/1B
AD-BASE

IT IS REQUIRED TO PRINT EACH NAME AT LEAST ONCE ON THIS PAGE

Released by:	To	Date	Purpose
	Andreas Brockbach	MAY 24 2004	begin Assay
	FREEZER #2	MAY 24 2004	STORAGE
FREEZER #2		MAY 26 2004	Continue Assay
	REFRIGERATOR 22	MAY 25 2004	Incubation
REFRIGERATOR 22		MAY 26 2004	Continue Assay
	Long term STORAGE	MAY 26 2004 2002-9-6 AVM AS SWH/T	SS. injection End Assay AS 05/26/04

UCLA LOGS & FORMS USED

Form approved SS. injection



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2535 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

CHAIN OF CUSTODY DOCUMENTATION

Code: BOTTLE OR ALIQUOT(S) A OR B

Aliquot: SCREEN (BATCH OR INDIVIDUAL) OR INDIVIDUAL CONFIRMATION

Organization: USADA UCLA folder no. USADA [REDACTED]

If batch, UCLA code numbers: _____

If individual, bottle no. [REDACTED] = UCLA code no. [REDACTED]

If multiple replicate aliquots indicate here N = _____

Circle assay: I B CF FF 2B TB ONTRAK DRU OTHER MS EPC Stryker
AC/BASE

IT IS REQUIRED TO PRINT EACH NAME AT LEAST ONCE ON THIS PAGE

Released by	To:	Date	Purpose
 Fareh Alhar Dehshad	 Andras Foldbach	MAY 24 2004	Begin Assay
	Water bath 2125	MAY 24 2004	Incubation
	Water bath 2125	MAY 25 2004	Continue Assay
	REFRIGERATOR 22	MAY 25 2004	Incubation
	REFRIGERATOR 22	MAY 26 2004	Continue Assay
	Long term STORAGE	MAY 26 2004	End Assay SI-5126104 05/26/04

UCLA OAL/USADA FORM 1250 (REV. 11/03)

Form approved:



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

DESCRIPTION OF METHODS

Detection of Recombinant Erythropoietic Proteins

SCREENING PROCEDURE

Urine samples that are tested for erythropoietic proteins are screened first. If the screening results suggests that a prohibited substance is present the samples are subjected to an A-confirmation procedure. In the UCLA laboratory the screen and the A-confirmation are performed by the same method. The main difference between the screen and the A-confirmation is that in the former several samples are electrophoresed on one gel next to each other, whereas in the latter the samples are separated by extra space.

The procedure described below detects recombinant human erythropoietin (rHuEPO) and darbepoetin. The analysis consists of four steps: sample preparation; isoelectric focusing (IEF), immuno-blotting, and visualization. The following description is intended to be a summary for the non-specialist.

Sample preparation

Preparing samples for IEF analysis involves two steps: 1) enzyme deactivation, 2) filtration, and concentration. The first step is to deactivate enzymes that could destroy the rHuEPO thereby foiling the analysis. It is done by adding protease inhibitors (chemicals) to the urine. The second step is to remove materials from the urine that are irrelevant and to concentrate the proteins. This is accomplished with specialized filters that retain molecules with high molecular weight, such as rHuEPO. The urine is placed in a cup that has the specialized filter. The cup is placed in a centrifuge which spins at high speed. This forces the low molecular weight material to pass through the filter. The rHuEPO and other proteins with similar molecular weights are retained (the retentate) on the filter. This step is repeated one more time. If rHuEPO is present in the urine, it will be in the final retentate. The final retentate is a liquid. A

small portion of it is used to estimate the concentration of rHuEPO by an immunoassay. Another portion is spotted on a gel (described below). The immunoassay provides an estimate of the amount of rHuEPO in the urine sample. This number is used to adjust the concentration of rHuEPO in the retentate to an optimal value. No adjustment is needed if the estimated concentration is low. If the estimate is high the retentate is diluted. Twenty microliters of the retentate are applied to the gel.

Isoelectric focusing

Glycoproteins such as darbepoetin and rHuEPO are molecules that carry positive and negative charges. The net charge of the glycoprotein is the algebraic sum of all the positive and negative charges. Each molecule has one net charge which is positive, negative, or neutral depending on the pH. The isoelectric point (pI) is a fundamental characteristic of proteins. It is the pH value at which the molecule is electrically neutral because the number of positive charges on the molecules are exactly balanced by the number of negative charges.

The gel is a jelly-like material that serves as the 'platform' for the electrophoresis. The gel is about 25 by 12 cm (length and width) and about 1 mm thick. The gel is prepared by mixing various chemicals in a flask and pouring the mixture into a cast. Just before the pour, a reagent is added that causes the materials to 'gel.' After the gel 'hardens' to the consistency of a flexible jelly-like material, the cast is removed and the gel is placed flat on the surface of the electrophoresis instrument.

Before the samples are put on the gel, a pH gradient must be set up. Electrodes are attached to the gel and connected to the electrophoresis unit. One of the electrodes attached to the gel is the anode and the other is the cathode. When the electrophoresis unit is turned on it sets up an electric circuit between the cathode (negative pole), and the anode (positive pole). When the current is applied the molecules used to set up the pH gradient migrate. The charged molecules migrate in the direction of the electrode bearing the opposite charge. Thus negatively charged or

more acidic molecules migrate toward the anode - the electrode with the positive charge and vice versa. A current is applied to the gel for half an hour to establish the pH gradient.

The current is turned off. The samples and standards that are to be electrophoresed are 'spotted' onto the gel by adding a small volume (20 microliters) of each sample to a piece of filter paper (10 by 5 mm) that has been placed on the gel. The pieces are placed 1 cm apart close to one edge of the gel. The gel can accommodate about 24 such pieces. In this way the surface of the gel is divided into 24 imaginary lanes.

Each sample or standard is spotted in one lane. The number of samples and standards that are processed on one gel is determined by the intent of the analysis. Typically we use one or two different standards, one or more control samples (content known), and several unknowns (content not known). Typically the standard is pure rHuEPO, or pure darbepoetin, or a mixture of the two. A typical control urine is a urine obtained from a subject to whom we have administered rHuEPO.

The current is turned back on which causes the glycoproteins that have been spotted on the gel to begin to move or migrate. The total number of hours is selected such that all the molecules have sufficient time to migrate or move to their isoelectric point. Once they reach the isoelectric point they remain stationary.

With this background, one can now interpret an electropherogram that shows rHuEPO (see figure 1). rHuEPO has 5 bands which are referred to as isoforms. An isoform is a subset of the rHuEPO molecules that has a defined pI. Not all the molecules of rHuEPO have exactly the same chemical structure and therefore the same pI. All the molecules in one band will have the same pI, that is they will become neutral at one and the same pH. The molecules of darbepoetin focus near the anode side of the gel because it has a low pI. In contrast, rHuEPO focuses closer to the cathode because it has a higher pI. This explains why rHuEPO and darbepoetin are separated on the gel by several centimeters of physical space.



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

Immuno-blotting

At this stage the bands are separated on the gel. The task now is to transfer them off the gel. This is accomplished by 'blotting'. Blotting is a procedure for transferring proteins from one surface to another. Two blotting procedures are used (blot 1 and blot 2). The first blot transfers the rHuEPO from the gel to membrane #1. Membrane #1 is incubated with antibodies against EPO proteins. The second blot transfers the antibodies from membrane #1 to membrane #2.


The gel is removed from the electrophoresis plate, washed with a buffer, and placed between two stacks of paper that have been soaked in a special blotting buffer. This stack (or sandwich) is placed into an instrument which is designed to apply an electrical charge across two plates. The stack is placed between the two plates and the current is turned on for 30 minutes. The instrument is referred to as the 'Blotting unit.' During blot #1 the rHuEPO 'travels' from the gel to membrane #1. Membrane #1 is a mirror image of the material that was on the gel.

Next membrane #1 is bathed in a solution of antibodies against EPO. These antibodies are very special because they specifically bind to endogenous human EPO, rHuEPO, and closely related proteins such as darbepoetin. The specific EPO/darbepoetin antibodies are obtained from mice that have been immunized and therefore make antibodies that react with darbepoetin and rHuEPO. Because we also use another antibody later in the assay, we sometimes refer to this antibody the 'first' or primary antibody.

In the next stage, referred to as the second blot, the primary antibody is transferred from the first membrane to a second membrane (membrane #2). The molecules of EPO and darbepoetin remain on the first membrane, but the antibodies against them (primary antibody) are transferred to the second membrane. The test detects the primary antibodies that mark the location of the darbepoetin and rHuEPO. The second blot is accomplished like the first blot. Then the second membrane is



UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



Incubated in a solution containing a second antibody that binds specifically to the primary antibody.

Visualization

The location of the second antibody on membrane #2 matches the location of EPO or darbepoetin on membrane #1. We need a way to visualize the second antibody. To this end, a marker protein which binds to the second antibody is used. Next a special substance is added that emits light when it comes in contact with the marker protein. The emitted light is then captured with a special digital camera. The final image (electropherogram) is used to evaluate the results.

CONFIRMATION PROCEDURE

The confirmation procedure is very similar to the screening procedure. The sample preparation, isoelectric focusing (IEF), immuno-blotting, and visualization are identical to the screening procedure. The only difference is the types and number of samples spotted on the confirmation gel. For the A-confirmation four lanes contain a standard (darbepoetin and rHuEPO), and there is one lane each for the Positive Quality Control, the Negative Quality Control, and the sample that is being confirmed.




Figure 1: Screen of the A sample

Electropherogram corresponding to screening a batch of urine samples for erythropoietic proteins. Each sample is contained in one lane. There are six lanes. The bands or isoforms corresponding to rHuEPO, darbepoetin (NESP), or endogenous EPO have the appearance of rungs of a ladder. There are two lanes (lanes 1 and 6, numbering from left to right) containing a standard consisting of a mixture of rHuEPO and darbepoetin, lane 2 contains a negative quality control urine sample, and lane 3 a positive quality control urine sample. Lane 4 is a sample from another athlete and lane 5 is sample # [REDACTED]



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

SAMPLE A SCREEN CONCLUSION

The two-band ratio is greater than 1.19 and the percent basic isoforms is greater than 80%, therefore an A-confirmation is carried out. (For the origin of the number 1.19 see publication: Breidbach A, Catlin DH, Green GA, Tregub I, Truong H, Gorzek J. Detection of rHuEPO in Urine by isoelectric focusing. Clin Chem 2003; 49:901-907)




Figure 2: Electropherogram corresponding to the A confirmation of sample # [REDACTED]

Each sample is contained in one lane. There are 7 lanes in the figure. The bands or isoforms corresponding to rHuEPO, darbepoetin (NESP), or endogenous EPO have the appearance of rungs of a ladder. There are four lanes containing a standard consisting of a mixture of rHuEPO and darbepoetin. These are lanes 1, 4, 5, and 7 with the numbering starting from the left-most lane. The negative QC sample is in lane 2 and the positive QC sample is in lane 3. The lane containing urine # [REDACTED] is lane 6.




Figure 3: Graphic representation of the band pattern - A confirmation

Densitometry has been used to express the darkness (density) of the bands in Figure 2 as numerical values. Plotting these values as a function of the band number results in graphs showing the relative peak height of the bands on the Y-axis and band number on the X-axis. The band labeled 0, which is marked by a vertical line, is defined as the band with the same isoelectric point (pI) as the most acidic band in the rHuEPO standard. Starting from band 0, the successive bands in the direction of the anode are considered 'acidic' and are labeled 1, 2, 3, 4,... Similarly the bands in the direction of the cathode are considered 'basic' and are labeled -1, -2, -3, and -4.

The shape of the graph for the negative QC is characteristic of a negative sample. The peak apex occurs at band 2 and the area under the curve is predominantly to the right of the vertical line. The shape of the graph for the positive QC is characteristic of a positive sample. The peak apex occurs at band -1 and the area under the curve is predominantly to the left of the vertical line.

For sample # [REDACTED] the peak apex is on the left of the vertical line like that of the positive QC and the area under the curve is predominantly on the left of the vertical line.



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

SAMPLE A CONFIRMATION CONCLUSIONS

1) The two-band ratio:

The A confirmation data indicate the presence of rHuEPO according to criteria published in a peer-reviewed journal by the UCLA Olympic Laboratory (Breibach A, Catlin DH, Green GA, Tregub I, Truong H, Gorzek J. Detection of rHuEPO in Urine by isoelectric focusing. Clin Chem 2003; 49:901-907.) In summary these criteria are:

- 1) the isoform bands that focus in the same area as the rHuEPO standard are darker than other isoform bands of the sample.
- 2) these bands have the same pI as the corresponding bands of the rHuEPO standard;
- 3) the sample contains the three bands that correspond to the pI of bands 0, -1, and -2 of the rHuEPO standard (see description for figure 3);
- 4) the 'Two-Band Ratio' is greater than 7. The upper 99% confidence interval for the two-band ratio is 1.19. (Note: the ratio is not a linear function (see Breidbach et.al.Clin Chem 2003; 49:901-907)

2) Percent basic isoforms:

In addition, the sample is positive for rHuEPO according to the 'percent basic isoforms' criteria. The percent basic isoforms is 90%. This method has been the subject of discussion in various expert committees hosted by the IOC and it has been described in an extended abstract by Pascual et. al. (Schaenzer et.al.(eds.) Recent Advances in Doping Analysis/ Proceedings of the Manfred Donike Workshop. 2002:135-144). In addition the method is discussed in a Court of Arbitration for Sport decision (CAS 2001/A/345 Roland Meier v/ Swiss Cycling) wherein a value greater than 80% was deemed to be positive



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01




Figure 4: Electropherogram corresponding to the B confirmation of sample # [REDACTED]

Each sample is contained in one lane. There are 7 lanes in the figure. The bands or isoforms corresponding to rHuEPO, darbepoetin (NESP), or endogenous EPO have the appearance of rungs of a ladder. There are four lanes containing a standard consisting of a mixture of rHuEPO and darbepoetin. These are lanes 1, 4, 5, and 7 with the numbering starting from the left-most lane. The negative QC sample is in lane 2 and the positive QC sample is in lane 3. The lane containing urine # [REDACTED] is lane 6.



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01




Figure 5 : Graphic representation of the band pattern - B confirmation

Densitometry has been used to express the darkness (density) of the bands in Figure 2 as numerical values. Plotting these values as a function of the band number results in graphs showing the relative peak height of the bands on the Y-axis and band number on the X-axis. The band labeled 0, which is marked by a vertical line, is defined as the band with the same isoelectric point (pI) as the most acidic band in the rHuEPO standard. Starting from band 0, the successive bands in the direction of the anode are considered 'acidic' and are labeled 1, 2, 3, 4,... Similarly the bands in the direction of the cathode are considered 'basic' and are labeled -1, -2, -3, and -4.

The shape of the graph for the negative QC is characteristic of a negative sample. The peak apex occurs at band 2 and the area under the curve is predominantly to the right of the vertical line. The shape of the graph for the positive QC is characteristic of a positive sample. The peak apex occurs at band -1 and the area under the curve is predominantly to the left of the vertical line.

For sample # [REDACTED] the peak apex is on the left of the vertical line like that of the positive QC and the area under the curve is predominantly on the left of the vertical line.



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

SAMPLE B CONFIRMATION CONCLUSIONS

1) The two-band ratio:

The A confirmation data indicate the presence of rHuEPO according to criteria published in a peer-reviewed journal by the UCLA Olympic Laboratory (Breidbach A, Catlin DH, Green GA, Tregub I, Truong H, Gorzek J. Detection of rHuEPO in Urine by isoelectric focusing. Clin Chem 2003; 49:901-907.) In summary these criteria are:

- 1) the isoform bands that focus in the same area as the rHuEPO standard are darker than other isoform bands of the sample.
- 2) these bands have the same pI as the corresponding bands of the rHuEPO standard;
- 3) the sample contains the three bands that correspond to the pI of bands 0, -1, and -2 of the rHuEPO standard (see description for figure 3);
- 4) the 'Two-Band Ratio' is greater than 6. The upper 99% confidence interval for the two-band ratio is 1.19. (Note: the ratio is not a linear function (see Breidbach et al. Clin Chem 2003; 49:901-907)

2) Percent basic isoforms:

In addition, sample # [REDACTED] is positive for rHuEPO according to the 'percent basic isoforms' criteria. The percent basic isoforms is 89%. This method has been the subject of discussion in various expert committees hosted by the IOC and it has been described in an extended abstract by Pascual et. al. (Schaenzer et.al.(eds.) Recent Advances in Doping Analysis/ Proceedings of the Manfred Donike Workshop. 2002:135-144). In addition the method is discussed in a Court of Arbitration for Sport decision (CAS 2001/A/345 Roland Meier v/ Swiss Cycling) wherein a value greater than 80% was deemed to be positive




IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077




ISO/IEC 17025
Biological Testing
Certificate: 1420-01

To verify that the isoforms are stable in the sample, we adjusted the pH to 5 with 0.1 mol/L HCl. We added Pepstatin (15 μmol/L to 500 μL) and Complete stock solution (1% of a 1 tablet/2mL solution) to 0.5 mL of sample. (Pepstatin and Complete are protease inhibitors.). After 15 min at room temperature 300mU rHuEPO were added and the



rHuEPO/NESP
USADA469853
Stability test
rHuEPO/NESP

Figure 6:
Stability
test - A
confirmati
on



rHuEPO/NESP
USADA469853
Stability test
rHuEPO/NESP

Figure 7:
Stability test -
B confirmation

urine was incubated at 37°C overnight. 20 μL of the treated urine sample were subjected to isoelectric focusing and immunoblotting. This study was carried out with sample # [redacted]. The rationale is if there is anything that causes a shift of the isoform pattern of the urine sample the well-known pattern of the added rHuEPO would change.

Figures 6 and 7 show the electropherograms of the above experiments. It can be seen that the isoform pattern of the added rHuEPO did not change.

Figure 8 shows an overlay of the isoform patterns of the screen, the A-, and the B-confirmation of sample # [REDACTED]. This further verifies that there is no change in the location of the isoforms between the three analyses.




Figure 8: Overlay of the graphic representations of the band pattern of screen, A-, and B-confirmation of sample # [REDACTED]

UNIVERSITY OF CALIFORNIA, LOS ANGELES

UCLA

BERKELEY • DAVIS • IRVINE • LOS ANGELES • RIVERSIDE • SAN DIEGO • SAN FRANCISCO

SANTA BARBARA • SANTA CRUZ



UCLA OLYMPIC ANALYTICAL LABORATORY
DEPARTMENT OF PHARMACOLOGY
UCLA SCHOOL OF MEDICINE
3123 GRANVILLE AVENUE
LOS ANGELES, CALIFORNIA 90025
PHONE (310) 825-2633
FAX (310) 206-9077

CONFIDENTIAL

June 21, 2004

Terry Madden
The United States Anti-Doping Agency
2550 Tenderfoot Hill St., Suite 200
Colorado Springs, CO 80906-7346

RE: Specimen number [REDACTED]

Dear Mr. Madden

Please find enclosed an addendum to the documentation package for the case identified above.

Please feel free to call if you have any questions.

Sincerely,

Don H. Catlin, M.D.
Director

**CONFIDENTIAL
DOCUMENTATION**

SAMPLE IDENTIFICATION:

Organization requesting test: USADA

Date of sample collection: Mar 16, 2004

Site ID: OOC



Substance identified: recombinant human Erythropoietin (rHuEPO)

TABLE OF CONTENTS

"A" SAMPLE CONFIRMATION DOCUMENTATION

Table 1 : Results of densitometry 3

"B" SAMPLE CONFIRMATION DOCUMENTATION

Table 2 : Results of densitometry 4



IOC Accredited

UCLA Olympic Analytical Laboratory
 UCLA School of Medicine
 2122 Granville Ave Los Angeles CA 90025
 Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
 Biological Testing
 Certificate: 1420-01

Band	Peak area
-4	37357
-3	232807
-2	463736
-1	553275
0	382344
1	99292
2	39495
3	24331
4	14888
5	10689

Table 1: Peak areas of the A confirmation of this sample

Densitometry has been used to convert the band densities (darkness) into areas (see documentation package page 27). From these areas the 'Two-band ratio' and the 'Percent basic isoforms' have been calculated. For the 'Two-band ratio' the sum of the areas of bands '-1' and '-2' is divided by the sum of the areas of bands '1' and '2'.

The value for the TBR is 7.3 $((553275+463736)/(99292+39495))$.

For the 'Percent basic isoforms' the sum of the areas of the basic bands ('0', '-1', '-2', etc.) is divided by the sum of the areas of all bands and then multiplied by 100.

The value for the 'Percent basic isoforms' is 89.8 %.



IOC Accredited

UCLA Olympic Analytical Laboratory
 UCLA School of Medicine
 2122 Granville Ave Los Angeles CA 90025
 Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
 Biological Testing
 Certificate: 1420-01

Band	Peak area
-4	22160
-3	204283
-2	510906
-1	665432
0	522979
1	125182
2	47359
3	29572
4	21253
5	12197

Table 2: Peak areas of the B confirmation of this sample

Densitometry has been used to convert the band densities (darkness) into areas (see documentation package page 27). From these areas the 'Two-band ratio' and the 'Percent basic isoforms' have been calculated. For the 'Two-band ratio' the sum of the areas of bands '-1' and '-2' is divided by the sum of the areas of bands '1' and '2'.

The value for the TBR is 6.8 $((665432+510906)/(125182+47359))$.

For the 'Percent basic isoforms' the sum of the areas of the basic bands ('0', '-1', '-2', etc.) is divided by the sum of the areas of all bands and then multiplied by 100.

The value for the 'Percent basic isoforms' is 89.1 %.



UCLA OLYMPIC ANALYTICAL LABORATORY
DEPARTMENT OF PHARMACOLOGY
UCLA SCHOOL OF MEDICINE
2121 GRANVILLE AVENUE
LOS ANGELES, CALIFORNIA 90024
PHONE (310) 825-2135
FAX (310) 206-9077

CONFIDENTIAL

November 12, 2004

Terry Madden
United States Anti-Doping Agency
2550 Tenderfoot Hill St., Suite 200
Colorado Springs, CO 80906

RE: Specimen number [REDACTED]

Dear Mr. Madden,

In response to your fax of today, please find enclosed an addendum to the documentation package for the case identified above.

In our publication (Breidbach A, Catlin DH, Green GA, Tregub I, Truong H, Gorzek J. Detection of rHuEPO in Urine by isoelectric focusing. Clin Chem 2003; 49:901-907) we state that a TBR value greater than 1.19 is associated with a risk of 1 in 100. Since that publication we have increased our rHuEPO-free control population to 685 samples, and based on that data, a TBR value of greater than 1.8 is associated with a risk of 1 in 100,000. That is, if 100,000 samples were analyzed we would expect to find one sample with a TBR greater than 1.8.

Please feel free to call if you have any questions

Sincerely,

Don H. Callin, M.D.
Director

**CONFIDENTIAL
DOCUMENTATION**

SAMPLE IDENTIFICATION:

Organization requesting test: USADA
Date of sample collection: Mar 16, 2004
Site ID: OOC

[REDACTED]
Substance identified: recombinant human Erythropoietin (rHuEPO)

TABLE OF CONTENTS

CRITERION FOR POSITIVITY	3
"A" SAMPLE CONFIRMATION DOCUMENTATION	
Figure 1 : Electropherogram with isoform patterns and band labeling according to WADA TD2004EPO draft of Oct 13, 2004	4
Figure 2 : Densitogram of this sample	5
"B" SAMPLE CONFIRMATION DOCUMENTATION	
Figure 3 : Electropherogram with isoform patterns and band labeling according to WADA TD2004EPO draft of Oct 13, 2004	6
Figure 4 : Densitogram of this sample	7



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01

CRITERION FOR POSITIVITY

As stated in the documentation package for this sample, dated Apr 15, 2004, our primary criterion for positivity is the "two-band ratio" (TBR). The concept of this ratio has undergone peer review and was published in 2003 (Breidbach A, Catlin DH, Green GA, Tregub I, Truong H, Gorzek J. Detection of rHuEPO in Urine by isoelectric focusing. Clin Chem 2003; 49:901-907).



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01




Figure 1: Electropherogram corresponding to the A confirmation of this sample.

Densitometry has been used to convert the band intensities (darkness) into numerical values (see documentation package page 27). Those numerical values (peak heights) are displayed next to the individual bands in parentheses. Also shown are the band IDs (1, 2, 3, 4,...) according to the World Anti-Doping Agency's (WADA) technical document TD2004EPO draft of Oct 13, 2004.

It can be seen that band 2 is the most intense band in the entire lane, that the second most intense band is band 3, and that there are three consecutive, acceptable bands (bands 1, 2, and 3) in the "basic" area. Therefore the identification criteria of WADA TD2004EPO draft of Oct 13, 2004 are satisfied.



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01




Figure 2 :Densitogram of lane USADA [redacted] in previous figure



IOC Accredited

UCLA Olympic Analytical Laboratory
UCLA School of Medicine
2122 Granville Ave Los Angeles CA 90025
Phone (310) 825-2635 FAX (310) 206-9077



ISO/IEC 17025
Biological Testing
Certificate: 1420-01




Figure 3: Electropherogram corresponding to the B confirmation of this sample

Densitometry has been used to convert the band intensities (darkness) into numerical values (see documentation package page 27). Those numerical values (peak heights) are displayed next to the individual bands in parentheses. Also shown are the band IDs (1, 2, 3, 4....) according to the World Anti-Doping Agency's (WADA) technical document TD2004EPO draft of Oct 13, 2004.

It can be seen that band 2 is the most intense band in the entire lane, that the second most intense band is band 1, and that there are three consecutive, acceptable bands (bands 1, 2, and 3) in the "basic" area. Therefore the identification criteria of WADA TD2004EPO draft of Oct 13, 2004 are satisfied.




Figure 4 : Densitogram of Jane USADA [REDACTED] in previous figure